





BCAP31 Antibody

Product Code	CSB-PA000973
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51572
Immunogen	Synthesized peptide derived from the Internal region of Human BAP31.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BCAP31; BAP31; DXS1357E; B-cell receptor-associated protein 31; BCR-associated protein 31; Bap31; 6C6-AG tumor-associated antigen; Protein CDM; p28
Immunogen Species	Homo sapiens (Human)
Target Names	BCAP31

Image



Western Blot analysis of 293T cells using BAP31 Polyclonal Antibody



Western Blot analysis of 293T AD293 22RV1
HELA HepG2 cells using BAP31 Polyclonal
Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.