





# BCAP31 Antibody

<b>Product Code</b>	CSB-PA000973
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P51572
<b>Immunogen</b>	Synthesized peptide derived from the Internal region of Human BAP31.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	BCAP31; BAP31; DXS1357E; B-cell receptor-associated protein 31; BCR-associated protein 31; Bap31; 6C6-AG tumor-associated antigen; Protein CDM; p28
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	BCAP31

## Image



Western Blot analysis of 293T cells using BAP31 Polyclonal Antibody



Western Blot analysis of 293T AD293 22RV1  
HELA HepG2 cells using BAP31 Polyclonal  
Antibody