





BCKDK Antibody

Product Code	CSB-PA000981
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14874
Immunogen	Synthesized peptide derived from the N-terminal region of Human BCKDK.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BCKDK; [3-methyl-2-oxobutanoate dehydrogenase [lipoamide]] kinase; mitochondrial; Branched-chain alpha-ketoacid dehydrogenase kinase; BCKD-kinase; BCKDHKIN
Immunogen Species	Homo sapiens (Human)
Target Names	BCKDK

Image



Western Blot analysis of 3T3 cells using BCKDK Polyclonal Antibody



Western Blot analysis of NIH-3T3 K562 Jurkat cells using BCKDK Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.