





BCL2 Antibody

Product Code	CSB-PA000982
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10415
Immunogen	Synthesized peptide derived from Human Bcl-2 around the non-phosphorylation site of S70.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BCL2; Apoptosis regulator Bcl-2
Immunogen Species	Homo sapiens (Human)
Target Names	BCL2

Image



Western Blot analysis of COS7 cells using Bcl-2 Polyclonal Antibody



Western Blot analysis of 293T cells using Bcl-2 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.