





FGFR2 Antibody

Product Code	CSB-PA000992
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21802
Immunogen	Synthesized peptide derived from the Internal region of Human Bek.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FGFR2; BEK; KGFR; KSAM; Fibroblast growth factor receptor 2; FGFR-2; K-sam; KGFR; Keratinocyte growth factor receptor; CD antigen CD332
Immunogen Species	Homo sapiens (Human)
Target Names	FGFR2

Image



Western Blot analysis of A549 cells using Bek Polyclonal Antibody



Western Blot analysis of HuvEc cells using Bek Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.