




BLK Antibody

Product Code	CSB-PA000997
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51451
Immunogen	Synthesized peptide derived from Human Blk around the non-phosphorylation site of Y501.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BLK; Tyrosine-protein kinase Blk; B lymphocyte kinase; p55-Blk
Immunogen Species	Homo sapiens (Human)
Target Names	BLK

Image



Western Blot analysis of 293T cells using Blk Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.