




BLNK Antibody

Product Code	CSB-PA000998
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WV28
Immunogen	Synthesized peptide derived from Human BLNK around the non-phosphorylation site of Y96.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BLNK; BASH; SLP65; B-cell linker protein; B-cell adapter containing a SH2 domain protein; B-cell adapter containing a Src homology 2 domain protein; Cytoplasmic adapter protein; Src homology 2 domain-containing leukocyte protein of 65 kDa;
Immunogen Species	Homo sapiens (Human)
Target Names	BLNK

Image



Western Blot analysis of 293 cells using BLNK Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.