





BMP3 Antibody

Product Code	CSB-PA001013
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12645
Immunogen	Synthesized peptide derived from the Internal region of Human BMP-3A.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BMP3; BMP3A; Bone morphogenetic protein 3; BMP-3; Bone morphogenetic protein 3A; BMP-3A; Osteogenin
Immunogen Species	Homo sapiens (Human)
Target Names	BMP3

Image



Western Blot analysis of HepG2 cells using BMP-3A Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using BMP-3A Polyclonal Antibody