





BMP3 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA001013 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P12645 |
| Immunogen | Synthesized peptide derived from the Internal region of Human BMP-3A. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | BMP3; BMP3A; Bone morphogenetic protein 3; BMP-3; Bone morphogenetic protein 3A; BMP-3A; Osteogenin |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | BMP3 |

Image



Western Blot analysis of HepG2 cells using BMP-3A Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using BMP-3A Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.