




BMX Antibody

Product Code	CSB-PA001034
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51813
Immunogen	Synthesized peptide derived from Human Bmx around the non-phosphorylation site of Y40.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BMX; Cytoplasmic tyrosine-protein kinase BMX; Bone marrow tyrosine kinase gene in chromosome X protein; Epithelial and endothelial tyrosine kinase; ETK; NTK38
Immunogen Species	Homo sapiens (Human)
Target Names	BMX

Image



Western Blot analysis of HepG2 cells using Bmx Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.