




BUD31 Antibody

Product Code	CSB-PA001072
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41223
Immunogen	Synthesized peptide derived from the N-terminal region of Human BUD31.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BUD31; EDG2; Protein BUD31 homolog; Protein EDG-2; Protein G10 homolog
Immunogen Species	Homo sapiens (Human)
Target Names	BUD31

Image



Western Blot analysis of Jurkat cells using BUD31 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.