




CA5A Antibody

Product Code	CSB-PA001099
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35218
Immunogen	Synthesized peptide derived from the Internal region of Human CA VA.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CA5A; CA5; Carbonic anhydrase 5A; mitochondrial; Carbonate dehydratase VA; Carbonic anhydrase VA; CA-VA
Immunogen Species	Homo sapiens (Human)
Target Names	CA5A

Image



Western Blot analysis of K562 cells using CA VA Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.