




CA13 Antibody

Product Code	CSB-PA001104
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N1Q1
Immunogen	Synthesized peptide derived from the Internal region of Human CA XIII.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CA13; Carbonic anhydrase 13; Carbonate dehydratase XIII; Carbonic anhydrase XIII; CA-XIII
Immunogen Species	Homo sapiens (Human)
Target Names	CA13

Image



Western Blot analysis of HepG2 cells using CA XIII Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.