




CAPN11 Antibody

Product Code	CSB-PA001185
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UMQ6
Immunogen	Synthesized peptide derived from the Internal region of Human Calpain 11.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CAPN11; Calpain-11; Calcium-activated neutral proteinase 11; CANP 11
Immunogen Species	Homo sapiens (Human)
Target Names	CAPN11

Image



Western Blot analysis of HuvEc cells using Calpain 11 Polyclonal Antibody