




CAPN11 Antibody

Product Code	CSB-PA001185
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UMQ6
Immunogen	Synthesized peptide derived from the Internal region of Human Calpain 11.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CAPN11; Calpain-11; Calcium-activated neutral proteinase 11; CANP 11
Immunogen Species	Homo sapiens (Human)
Target Names	CAPN11

Image



Western Blot analysis of HuvEc cells using Calpain 11 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.