




ANKRD1 Antibody

Product Code	CSB-PA001217
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15327
Immunogen	Synthesized peptide derived from the C-terminal region of Human CARP.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ANKRD1; C193; CARP; HA1A2; Ankyrin repeat domain-containing protein 1; Cardiac ankyrin repeat protein; Cytokine-inducible gene C-193 protein; Cytokine-inducible nuclear protein
Immunogen Species	Homo sapiens (Human)
Target Names	ANKRD1

Image



Western Blot analysis of COLO205 cells using CARP Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.