





CASP9 Antibody

Product Code	CSB-PA001235
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55211
Immunogen	Synthesized peptide derived from Human Caspase-9 around the non-phosphorylation site of T125.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CASP9; MCH6; Caspase-9; CASP-9; Apoptotic protease Mch-6; Apoptotic protease-activating factor 3; APAF-3; ICE-like apoptotic protease 6; ICE-LAP6
Immunogen Species	Homo sapiens (Human)
Target Names	CASP9

Image



Western Blot analysis of 3T3 cells using Caspase-9 Polyclonal Antibody



Western Blot analysis of L929 rat-muscle mouse-brain rat-kidney cells using Caspase-9 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.