




CTNNA1/CTNNA2 Antibody

Product Code	CSB-PA001267
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35221/P26232
Immunogen	Synthesized peptide derived from the C-terminal region of Human Catenin-α E/N.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CTNNA1; Catenin alpha-1; Alpha E-catenin; Cadherin-associated protein; Renal carcinoma antigen NY-REN-13; CTNNA2; CAPR; Catenin alpha-2; Alpha N-catenin; Alpha-catenin-related protein
Immunogen Species	Homo sapiens (Human)
Target Names	CTNNA1/CTNNA2

Image



Western Blot analysis of HeLa cells using Catenin-α E/N Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.