





CREBBP Antibody

Product Code	CSB-PA001359
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92793
Immunogen	Synthesized peptide derived from the C-terminal region of Human CBP.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CREBBP; CBP; CREB-binding protein
Immunogen Species	Homo sapiens (Human)
Target Names	CREBBP

Image



Western Blot analysis of HeLa cells using CBP Polyclonal Antibody



Western Blot analysis of 293T cells using CBP Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.