




ACKR4 Antibody

Product Code	CSB-PA001407
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPB9
Immunogen	Synthesized peptide derived from the N-terminal region of Human CCRL1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CCRL1; CCBP2; CCR11; VSHK1; C-C chemokine receptor type 11; C-C CKR-11; CC-CKR-11; CCR-11; CC chemokine receptor-like 1; CCRL1; CCX CKR
Immunogen Species	Homo sapiens (Human)
Target Names	ACKR4

Image



Western Blot analysis of 293 cells using CCRL1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.