




SIRPB1 Antibody

Product Code	CSB-PA001428
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00241
Immunogen	Synthesized peptide derived from the Internal region of Human CD172b.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SIRPB1; Signal-regulatory protein beta-1; SIRP-beta-1; CD172 antigen-like family member B; CD antigen CD172b
Immunogen Species	Homo sapiens (Human)
Target Names	SIRPB1

Image



Western Blot analysis of HuvEc cells using CD172b Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.