





TLR4 Antibody

Product Code	CSB-PA001434
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00206
Immunogen	Synthesized peptide derived from the Internal region of Human CD284.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TLR4; Toll-like receptor 4; hToll; CD antigen CD284
Immunogen Species	Homo sapiens (Human)
Target Names	TLR4

Image



Western Blot analysis of HeLa cells using CD284 Polyclonal Antibody



Western Blot analysis of KB HELA cells using CD284 Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.