





CD4 Antibody

Product Code	CSB-PA001453
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01730
Immunogen	Synthesized peptide derived from Human CD4 around the non-phosphorylation site of S433.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD4; T-cell surface glycoprotein CD4; T-cell surface antigen T4/Leu-3; CD antigen CD4
Immunogen Species	Homo sapiens (Human)
Target Names	CD4

Image



Western Blot analysis of COLO205 cells using CD4 Polyclonal Antibody



Western Blot analysis of HepG2 cells using CD4 Polyclonal Antibody