





CD44 Antibody

Product Code	CSB-PA001460
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16070
Immunogen	Synthesized peptide derived from Human CD44 around the non-phosphorylation site of S706.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD44; LHR; MDU2; MDU3; MIC4; CD44 antigen; CDw44; Epican; Extracellular matrix receptor III; ECMR-III; GP90 lymphocyte homing/adhesion receptor; HUTCH-I; Heparan sulfate proteoglycan; Hermes antigen; Hyaluronate receptor; Phagocytic glycopr
Immunogen Species	Homo sapiens (Human)
Target Names	CD44

Image



Western Blot analysis of hela cells using CD44 Polyclonal Antibody



Western Blot analysis of 293T HeLa 823 cells using CD44 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.