




CD5 Antibody

Product Code	CSB-PA001470
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06127
Immunogen	Synthesized peptide derived from Human CD5 around the non-phosphorylation site of Y453.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD5; LEU1; T-cell surface glycoprotein CD5; Lymphocyte antigen T1/Leu-1; CD antigen CD5
Immunogen Species	Homo sapiens (Human)
Target Names	CD5

Image



Western Blot analysis of 293 cells using CD5 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.