




CD58 Antibody

Product Code	CSB-PA001472
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19256
Immunogen	Synthesized peptide derived from the Internal region of Human CD58.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CD58; LFA3; Lymphocyte function-associated antigen 3; Ag3; Surface glycoprotein LFA-3; CD antigen CD58
Immunogen Species	Homo sapiens (Human)
Target Names	CD58

Image



Western Blot analysis of HuvEc cells using CD58 Polyclonal Antibody