




TFRC Antibody

Product Code	CSB-PA001477
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02786
Immunogen	Synthesized peptide derived from Human CD71 around the non-phosphorylation site of S24.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TFRC; Transferrin receptor protein 1; TR; TfR; TfR1; Trfr; T9; p90; CD antigen CD71
Immunogen Species	Homo sapiens (Human)
Target Names	TFRC

Image



Western Blot analysis of 293 cells using CD71 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.