




C5AR1 Antibody

Product Code	CSB-PA001479
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21730
Immunogen	Synthesized peptide derived from Human CD88 around the non-phosphorylation site of S338.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	C5AR1; C5AR; C5R1; C5a anaphylatoxin chemotactic receptor; C5a-R; C5aR; CD antigen CD88
Immunogen Species	Homo sapiens (Human)
Target Names	C5AR1

Image



Western Blot analysis of HeLa cells using CD88 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.