





FAS Antibody

Product Code	CSB-PA001488
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25445
Immunogen	Synthesized peptide derived from the C-terminal region of Human CD95.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FAS; APT1; FAS1; TNFRSF6; Tumor necrosis factor receptor superfamily member 6; Apo-1 antigen; Apoptosis-mediating surface antigen FAS; FASLG receptor; CD antigen CD95
Immunogen Species	Homo sapiens (Human)
Target Names	FAS

Image



Western Blot analysis of 823 cells using CD95 Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using CD95 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.