





MCM5 Antibody

Product Code	CSB-PA001513
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33992
Immunogen	Synthesized peptide derived from the N-terminal region of Human CDC46.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MCM5; CDC46; DNA replication licensing factor MCM5; CDC46 homolog; P1-CDC46
Immunogen Species	Homo sapiens (Human)
Target Names	MCM5

Image



Western Blot analysis of HepG2 cells using CDC46 Polyclonal Antibody



Western Blot analysis of NIH-3T3 Jurkat cells
using CDC46 Polyclonal Antibody