




CNTROB Antibody

Product Code	CSB-PA001554
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N137
Immunogen	Synthesized peptide derived from the Internal region of Human Centrobin.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CNTROB; LIP8; PP1221; Centrobin; Centrosomal BRCA2-interacting protein; LYST-interacting protein 8
Immunogen Species	Homo sapiens (Human)
Target Names	CNTROB

Image



Western Blot analysis of HuvEc cells using Centrobin Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.