





CEP170 Antibody

Product Code	CSB-PA001570
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5SW79
Immunogen	Synthesized peptide derived from the Internal region of Human CEP170.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CEP170; FAM68A; KAB; KIAA0470; Centrosomal protein of 170 kDa; Cep170; KARP-1-binding protein; KARP1-binding protein
Immunogen Species	Homo sapiens (Human)
Target Names	CEP170

Image



Western Blot analysis of AD293 cells using CEP170 Polyclonal Antibody



Western Blot analysis of MCF7 PC-3 HeLa cells using CEP170 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.