




CEP55 Antibody

Product Code	CSB-PA001584
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53EZ4
Immunogen	Synthesized peptide derived from the Internal region of Human CEP55.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CEP55; C10orf3; URCC6; Centrosomal protein of 55 kDa; Cep55; Up-regulated in colon cancer 6
Immunogen Species	Homo sapiens (Human)
Target Names	CEP55

Image



Western Blot analysis of HuvEc cells using CEP55 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.