





DDIT3 Antibody

Product Code	CSB-PA001640
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35638
Immunogen	Synthesized peptide derived from Human CHOP around the non-phosphorylation site of S30.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DDIT3; CHOP; CHOP10; GADD153; DNA damage-inducible transcript 3 protein; DDIT-3; C/EBP-homologous protein; CHOP; C/EBP-homologous protein 10; CHOP-10; Growth arrest and DNA damage-inducible protein GADD153
Immunogen Species	Homo sapiens (Human)
Target Names	DDIT3

Image



Western Blot analysis of VEC cells using CHOP Polyclonal Antibody



Western Blot analysis of HepG2 cells using CHOP Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.