




CLDN19 Antibody

Product Code	CSB-PA001657
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6F1
Immunogen	Synthesized peptide derived from the Internal region of Human Claudin-19.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CLDN19; Claudin-19
Immunogen Species	Homo sapiens (Human)
Target Names	CLDN19

Image



Western Blot analysis of HepG2 cells using Claudin-19 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.