




CLCN7 Antibody

Product Code	CSB-PA001686
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51798
Immunogen	Synthesized peptide derived from the N-terminal region of Human CLC-7.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CLCN7; H(+)/Cl(-) exchange transporter 7; Chloride channel 7 alpha subunit; Chloride channel protein 7; CIC-7
Immunogen Species	Homo sapiens (Human)
Target Names	CLCN7

Image



Western Blot analysis of LOVO cells using CLC-7 Polyclonal Antibody