





CLK2 Antibody

Product Code	CSB-PA001718
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49760
Immunogen	Synthesized peptide derived from the N-terminal region of Human CLK2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CLK2; Dual specificity protein kinase CLK2; CDC-like kinase 2
Immunogen Species	Homo sapiens (Human)
Target Names	CLK2

Image



Western Blot analysis of A549 cells using CLK2 Polyclonal Antibody



Western Blot analysis of HepG2 cells using CLK2 Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.