




CLK4 Antibody

Product Code	CSB-PA001719
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HAZ1
Immunogen	Synthesized peptide derived from the Internal region of Human CLK4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CLK4; Dual specificity protein kinase CLK4; CDC-like kinase 4
Immunogen Species	Homo sapiens (Human)
Target Names	CLK4

Image



Western Blot analysis of HT29 cells using CLK4 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.