





MYC Antibody

Product Code	CSB-PA001729
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01106
Immunogen	Synthesized peptide derived from Human c-Myc around the non-phosphorylation site of T358.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MYC; BHLHE39; Myc proto-oncogene protein; Class E basic helix-loop-helix protein 39; bHLHe39; Proto-oncogene c-Myc; Transcription factor p64
Immunogen Species	Homo sapiens (Human)
Target Names	MYC

Image



Western Blot analysis of HuvEc cells using c-Myc Polyclonal Antibody



Western Blot analysis of various cells using c-Myc Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.