




CNOT7 Antibody

Product Code	CSB-PA001734
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIV1
Immunogen	Synthesized peptide derived from the N-terminal region of Human CNOT7.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CNOT7; CAF1; CCR4-NOT transcription complex subunit 7; BTG1-binding factor 1; CCR4-associated factor 1; CAF-1
Immunogen Species	Homo sapiens (Human)
Target Names	CNOT7

Image



Western Blot analysis of HuvEc cells using CNOT7 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.