




GJB7 Antibody

Product Code	CSB-PA001754
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PEY0
Immunogen	Synthesized peptide derived from the Internal region of Human Connexin 25.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GJB7; CX25; Gap junction beta-7 protein; Connexin-25; Cx25
Immunogen Species	Homo sapiens (Human)
Target Names	GJB7

Image



Western Blot analysis of HT29 cells using Connexin 25 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.