




COPZ1 Antibody

Product Code	CSB-PA001759
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61923
Immunogen	Synthesized peptide derived from the N-terminal region of Human COP ζ1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	COPZ1; COPZ; CGI-120; HSPC181; Coatomer subunit zeta-1; Zeta-1-coat protein; Zeta-1 COP
Immunogen Species	Homo sapiens (Human)
Target Names	COPZ1

Image



Western Blot analysis of COS7 cells using COP ζ1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.