




COX17 Antibody

Product Code	CSB-PA001762
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14061
Immunogen	Synthesized peptide derived from the N-terminal region of Human COX17.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	COX17; Cytochrome c oxidase copper chaperone
Immunogen Species	Homo sapiens (Human)
Target Names	COX17

Image



Western Blot analysis of 293 cells using COX17 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.