





CSK Antibody

Product Code	CSB-PA001798
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41240
Immunogen	Synthesized peptide derived from Human Csk around the non-phosphorylation site of S364.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CSK; Tyrosine-protein kinase CSK; C-Src kinase; Protein-tyrosine kinase CYL
Immunogen Species	Homo sapiens (Human)
Target Names	CSK

Image



Western Blot analysis of 22RV1 cells using Csk Polyclonal Antibody



Western Blot analysis of Jurkat cells using Csk Polyclonal Antibody