





COPS3 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA001799 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UNS2 |
| Immunogen | Synthesized peptide derived from the Internal region of Human CSN3. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | COPS3; CSN3; COP9 signalosome complex subunit 3; SGN3; Signalosome subunit 3; JAB1-containing signalosome subunit 3 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | COPS3 |

Image



Western Blot analysis of 3T3 cells using CSN3 Polyclonal Antibody



Western Blot analysis of NIH-3T3 Jurkat RAT-intestine cells using CSN3 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.