




CUL3 Antibody

Product Code	CSB-PA001841
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13618
Immunogen	Synthesized peptide derived from the N-terminal region of Human CUL-3.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CUL3; KIAA0617; Cullin-3; CUL-3
Immunogen Species	Homo sapiens (Human)
Target Names	CUL3

Image



Western Blot analysis of Jurkat cells using CUL-3 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.