





CXCR3 Antibody

Product Code	CSB-PA001843
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49682
Immunogen	Synthesized peptide derived from the Internal region of Human CXCR-3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CXCR3; GPR9; C-X-C chemokine receptor type 3; CXCR-3; CKR-L2; G protein-coupled receptor 9; Interferon-inducible protein 10 receptor; IP-10 receptor; CD antigen CD183
Immunogen Species	Homo sapiens (Human)
Target Names	CXCR3

Image



Western Blot analysis of K562 cells using CXCR-3 Polyclonal Antibody



Western Blot analysis of various cells using CXCR-3 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.