





CCNB1 Antibody

Product Code	CSB-PA001863
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14635
Immunogen	Synthesized peptide derived from Human Cyclin B1 around the non-phosphorylation site of S126.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CCNB1; CCNB; G2/mitotic-specific cyclin-B1
Immunogen Species	Homo sapiens (Human)
Target Names	CCNB1

Image



Western Blot analysis of PC-12 cells using Cyclin B1 Polyclonal Antibody



Western Blot analysis of various cells using
Cyclin B1 Polyclonal Antibody