




CCNH Antibody

Product Code	CSB-PA001888
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51946
Immunogen	Synthesized peptide derived from Human Cyclin H around the non-phosphorylation site of T315.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CCNH; Cyclin-H; MO15-associated protein; p34; p37
Immunogen Species	Homo sapiens (Human)
Target Names	CCNH

Image



Western Blot analysis of HuvEc cells using Cyclin H Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.