





CYP1B1 Antibody

Product Code	CSB-PA001931
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16678
Immunogen	Synthesized peptide derived from the Internal region of Human CYP1B1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CYP1B1; Cytochrome P450 1B1; CYP1B1
Immunogen Species	Homo sapiens (Human)
Target Names	CYP1B1

Image



Western Blot analysis of Jurkat cells using CYP1B1 Polyclonal Antibody



Western Blot analysis of Jurkat cells using CYP1B1 Polyclonal Antibody

